

Accepted Manuscript

Title: High-performance liquid chromatographic determination of 2-aminoethylphosphonic acid and 2-amino-3-phosphonopropionic acid in seawater matrix using precolumn fluorescence derivatization with *o*-phthalaldehyde-ethanethiol



Authors: Shu Wang, Yao Deng, Baomin Liu, Xiaolin Li, Xin Lin, Dongxing Yuan, Jian Ma

PII: S0021-9673(18)31001-X
DOI: <https://doi.org/10.1016/j.chroma.2018.08.016>
Reference: CHROMA 359612

To appear in: *Journal of Chromatography A*

Received date: 12-4-2018
Revised date: 21-7-2018
Accepted date: 6-8-2018

Please cite this article as: Wang S, Deng Y, Liu B, Li X, Lin X, Yuan D, Ma J, High-performance liquid chromatographic determination of 2-aminoethylphosphonic acid and 2-amino-3-phosphonopropionic acid in seawater matrix using precolumn fluorescence derivatization with *o*-phthalaldehyde-ethanethiol, *Journal of Chromatography A* (2018), <https://doi.org/10.1016/j.chroma.2018.08.016>

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.

High-performance liquid chromatographic determination of 2-aminoethylphosphonic acid and 2-amino-3-phosphonopropionic acid in seawater matrix using precolumn fluorescence derivatization with *o*-phthalaldehyde-ethanethiol

Shu Wang^a, Yao Deng^a, Baomin Liu^b, Xiaolin Li^b, Xin Lin^b, Dongxing Yuan^a, Jian Ma^{a*}

^a State Key Laboratory of Marine Environmental Science, College of the Environment and Ecology, Xiamen University, 361102, Xiamen, Fujian, China

^b State Key Laboratory of Marine Environmental Science, College of the Ocean and Earth Sciences, Xiamen University, 361102, Xiamen, Fujian, China

*: Corresponding author: jma@xmu.edu.cn (Dr. J. Ma) Tel./Fax: 0086-592-2186916

Highlights

- Simple HPLC-FLD analysis 2-AEP and 2-AP3 after one-step OPA-ET derivatization.
- One of few methods for the determination of 2-AEP and 2-AP3 in seawater matrix.
- Careful evaluation of salinity effect on derivative efficiency and detection.
- Application to analysis of seawater algal culture media samples.

Download English Version:

<https://daneshyari.com/en/article/10154467>

Download Persian Version:

<https://daneshyari.com/article/10154467>

[Daneshyari.com](https://daneshyari.com)