



Supercritical fluid chromatographic resolution of water soluble isomeric carboxyl/amine terminated peptides facilitated via mobile phase water and ion pair formation

M.A. Patel^a, F. Riley^a, M. Ashraf-Khorassani^b, L.T. Taylor^{b,*}

^a Global Research & Development, Pfizer, Inc., Groton, CT 06340, USA

^b Department of Chemistry, Virginia Tech, Blacksburg, VA 24061-0212, USA

ARTICLE INFO

Article history:

Received 8 December 2011

Accepted 9 February 2012

Available online 17 February 2012

Keywords:

Supercritical fluid chromatography

Peptides

Isomers

Mass spectrometry

Evaporative light scattering

ABSTRACT

Both analytical scale and preparative scale packed column supercritical fluid chromatography (SFC) have found widespread applicability for chiral separations of multiple polar pharmaceutical candidates. However, SFC is rapidly becoming an achiral technique. More specifically, ion pair SFC is finding greater utility for separation of ionic analytes such as amine salts and organic sulfonates. The key to this success is, in part, the incorporation of additives such as trifluoroacetic acid and ammonium acetate into the mobile phase in association with a wide variety of both bonded silica stationary phases and high purity bare silica. Ion pairing SFC coupled with evaporative light scattering detection and mass spectrometric detection is presented here for the separation of water soluble, uncapped, isomeric peptide pairs that differ in amino acid arrangement. The separation is best achieved on either diol-bonded silica or bare silica with 1–5% (w/w) water as a significant ingredient in the mobile phase. Nitrogenous stationary phases such as 2-ethylpyridine, which had been very successful for the separation of capped peptides failed to yield the desired separation regardless of the mobile phase composition. A HILIC type retention mechanism is postulated for the separation of both isomeric uncapped peptide pairs.

© 2012 Elsevier B.V. All rights reserved.

1. Introduction

The use of mobile phase additives has been quite successful for separation of polar analytes via packed column supercritical fluid chromatography (pcSFC) [1,2]. The additives are generally added directly to the modifier and then pumped as a single binary fluid thereby augmenting the solvating power and flow of pressurized carbon dioxide. Additives are thought to perform at least five different functions [3]. These include: (a) cover active sites on the solid support, (b) change the polarity of the stationary phase, (c) suppress ionization of the analyte, (d) enhance ion pair formation with the analyte, and (e) raise the solvating power of the mobile phase. Of these five functions, modification of the stationary phase may be the most critical function. Maximum adsorption of CO₂ has been found to occur near its critical temperature and critical pressure [4]. Parcher et al. have suggested that a thick film of CO₂ may be coating the stationary phase and acting as a part of it. Modifiers have been shown to strongly adsorb to the stationary phase [5], and additives such as organic acids and bases, inorganic salts, and water are expected to do likewise.

Low to moderately polar stationary phases such as octyl and cyanopropyl were observed to adsorb only small amounts (0.4–0.6% of a monolayer) of modifier and additive. Under the same conditions, sulfonic acid and diol columns adsorbed much larger amounts of modifier/additive creating surface coverage up to 21% of a monolayer [6]. Basic, acidic, and salt additives have proven most popular, but as this report will attest, neutral water offers some significant advantages. Such behavior is not unique to SFC since it has been theorized and now demonstrated that a water-rich layer is associated with the silica surface in hydrophilic interaction chromatography (HILIC) [7].

Applications of ion-pairing SFC have been demonstrated in the past [8,9]. Pinkston et al. brought new thinking to this area when they showed that low levels of volatile ammonium salts as mobile phase additives allowed the elution of polar and even ionic organic materials from a Deltabond cyanopropyl-bonded silica column [10]. However, it was not clear from this study whether the salt ion paired with the analyte or if the salt interacted with the stationary phase. In a more recent study, the effect of various ammonium salts on the elution of *anionic* compounds such as the sodium salts of 4-dodecylbenzene sulfonate and 4-octylbenzene sulfonate has been shown [11]. The successful elution of *cationic* species (i.e. NR₄⁺X[−]), from a Deltabond cyanopropyl-bonded silica column was achieved by adding sodium alkylsulfonate to the methanol-modified carbon

* Corresponding author. Tel.: +1 540 231 6680; fax: +1 540 231 3255.

E-mail address: ltaylor@vt.edu (L.T. Taylor).

dioxide-based mobile phase. An ion pairing interaction between the positively charged amine cation and the anionic alkylsulfonate ion was proposed [12].

In some cases, conversion of a neutral analyte into an ion actually enhances SFC separation by favoring an ion-pairing mechanism. Stringham has reported that incorporation of a strong acid in both the sample diluent and mobile phase modifier gives a dramatic improvement in the separation of chiral basic compounds [13].

The re-discovery that the use of water as a neutral additive can lead to striking separations with bare silica has given new emphasis to the area of packed column SFC. Subcritical fluid chromatography of monosaccharides and polyols using bare silica and trimethylsilyl columns with a mobile phase composed of CO₂–methanol–water was reported nearly 15 years ago [14]. By adjusting the column temperature to 60 °C and the flow rate to 5 mL/min, a complete separation of eight underivatized monosaccharides was obtained in less than 10 min. Higher efficiencies and resolution were obtained when both water and triethylamine were added to methanol as the polar modifier.

Water, of course, is impractical as a direct modifier due to its very low solubility (~0.1%) in supercritical CO₂, but water blended as a tertiary component into alcohol modified supercritical CO₂ can assist in packed column SFC separations of highly polar analytes and thereby expand the capabilities of the technique. In this regard, 2-propanol and methanol have been compared as alcohol/water modifiers of CO₂ [15]. More recently, one study has described the packed column SFC separation of nucleobases on a variety of stationary phases with methanol/water as the modifier/additive [16]. In another study, numerous polar analytes and drug-like molecules have been shown to yield sharp peaks when eluted from bare silica with a mixture of methanol–water–CO₂ [17]. The removal of water from the modifier resulted in a decrease in selectivity and an increase in retention factor, however, the sharp peak shape was still achieved when water was removed from the mobile phase. The excellent chromatography could be restored via incorporation of water into the mobile phase. Another interesting feature of water as an additive is its compatibility with preparative SFC.

Successful pcSFC of polypeptides has employed a variety of mobile phase additives [18]. Blackwell and Stringham eluted peptides such as bradykinine and oxytocin from a divinylbenzene polymeric (i.e. 500 Å) column using heptadecafluorooctanesulfonic acid as an additive in methanol with CO₂ [19]. The separation of polypeptide salts containing a variety of acidic and basic residues was achieved via pcSFC–MS when trifluoroacetic acid (TFA) was used as the additive in a CO₂/methanol mobile phase on a 2-ethylpyridine–silica column [20].

With the discovery that ionic peptides of moderately high molecular mass can be eluted via pcSFC by incorporation of selected mobile phase additives, our interest has turned to the separation of polypeptide mixtures. A study was recently reported wherein the separation of two pairs of linear, water soluble, 12mer peptides of identical molecular mass (~1200 Da), composition, and charge that differ only in amino acid sequence was reported [21]. TFA as the additive in methanol proved to be the most successful modifier for elution of the two isomeric peptide pairs. A variety of silica-based basic stationary phases such as 3-aminopropyl and 2-ethylpyridine successfully eluted the isomeric peptide pairs. The decapeptide pairs in this case were end-capped wherein the C-terminus was an amide group and the N-terminus was an acetyl group. The lysine side-group was protonated but not capped.

We would now like to report the separation of water soluble, decapeptide pairs of identical mass, composition, and charge that differ only in amino acid sequence. In this study, the C and N terminus nor the side lysine group have been capped. The side-group however was again protonated. As in our previous

work referenced above, TFA in conjunction with methanol as the modifier proved to be successful for the elution of each of the carboxyl-terminated/amine-terminated peptide pairs via IPSFC. What is unique in this study is that the resolution of the isomeric unprotected peptide pairs was uniquely achieved when water was incorporated into the gradient mobile phase and either bare silica or diol-bonded silica was used as the stationary phase [22].

2. Experimental

2.1. Peptides and sample preparation

Uncapped peptides containing both an amine and carboxyl terminus were synthesized by CPC Scientific Inc. (San Jose, CA) for this research. The linear isomeric peptides exhibited high water solubility (e.g. 1 mg/mL). One pair of peptides (CI-05-00457 and CI-05-00444) exhibited a molecular mass of 1125.4 Da and the other pair of isomeric peptides (CI-05-00447 and CI-05-00448) had a molecular mass of 1173.5 Da. Each peptide pair had the same composition and exhibited minimal variation in the amino acid sequence. The former pair (1125.4) differed only in the placement of one of the glycines (G) and the lone valine (V) in the chain. The latter peptide pair (1173.5) differed only in the placement of one of the glycines (G) and the lone phenylalanine (F). Each peptide chain contained two lysine units with one nitrogen atom in the linear chain and the second nitrogen appearing as a dangling amino function that was no doubt positively charged in the acidic CO₂-containing mobile phase. The peptides supplied by CPC Scientific Inc. were in fact stated to be salts of trifluoroacetate, which served as the negative counterion. The amino acid sequence along with the quoted molecular mass and injected sample solution concentration for each of the amine/carboxyl terminated, cationic peptide chains are given in Table 1. Samples for chromatographic injections on the SFC–MS were prepared by pipetting a 50:50 mixture of each peptide pair, for example 0.5 mL of CI-05-00444 + 0.5 mL of CI-05-00457. Samples for chromatographic injections on the SFC–ELSD were prepared in a similar fashion although the total volumes were smaller.

2.2. Instrumentation and chromatographic conditions

All experiments performed at Virginia Tech were done using a Waters Corp. (Milford, MA) SFC equipped with a high pressure pump, evaporative light scattering detector (Alltech ELSD), auto-sampler, oven heater set to 60 °C, and back pressure regulator set at 140 bar. All SFC analyses were performed via gradient elution. All experiments performed at Pfizer were done using a Waters Analytical method station II with a ZQ mass spectrometer (Milford, MA). The back pressure regulator was set at 120 bar and the oven temperature was set at 60 °C. The percent additive in the methanol modifier is given in the figure captions. The flow rate was 2 mL/min and the linear gradient schedule for this study is given in Table 2. For ELSD measurements, the outlet flow was split post-back pressure restrictor (BPR) such that only 1/4 of the flow passed through the ELSD. The ELSD drift tube temperature was set to 105 °C and the nitrogen flow was set to approximately 3.5 L/min. The split ratio was adjusted with a fused silica restrictor. Sample injection volume was 5 µL and the concentration of the injected solution was approximately 0.5 µg/µL/peptide isomer. In addition to bare silica, a number of polar and non-polar stationary phases were examined as candidates for both efficient elution of each peptide and near baseline resolution of the isomeric carboxyl/amine terminated peptide salts under supercritical conditions. A list of all columns extensively investigated in this study can be found in Table 3.

Download English Version:

<https://daneshyari.com/en/article/1206370>

Download Persian Version:

<https://daneshyari.com/article/1206370>

[Daneshyari.com](https://daneshyari.com)