



Development and validation of an assay method for the determination of trifluoroacetic acid in a cyclosporin-like drug

Mark Powell^{a,*}, David Humphreys^a, Adam Nagle^a, Karen Polowy^{b,1}, Michael Scannell^{b,2}

^a Quay Pharmaceuticals Limited, 28 Parkway, Deeside, Flintshire CH5 2NS, United Kingdom

^b Scynexis Incorporated, 3501 C Tricenter Blvd., Durham, NC 27713, United States

ARTICLE INFO

Article history:

Received 16 February 2011

Received in revised form 28 March 2011

Accepted 11 April 2011

Available online 6 May 2011

Keywords:

Ion chromatography

Trifluoroacetic acid

Cyclosporin

ABSTRACT

An improved method for the assay of trifluoroacetic acid (TFA) in a cyclosporin-like drug substance is presented, based on ion chromatography with suppressed conductivity detection. Column fouling by the drug molecule is avoided by use of a sample preparation method in which the drug substance is precipitated at alkaline pH whilst the TFA remains in solution. The new method requires a smaller sample mass than a previous method based on headspace-GC-FID whilst achieving an improvement in sensitivity. During validation, the method's performance was found to be consistent with usual acceptance criteria, and the method was found to be robust in routine use.

© 2011 Elsevier B.V. All rights reserved.

1. Introduction

SCY-502635 (Fig. 1) is a novel non-immunosuppressive cyclosporin-based analogue under investigation for the treatment of chronic human hepatitis C infection (HCV; [1,2]). It is similar in structure to cyclosporin A; the differences arise through modification of the residues at the 3 and 4 positions of the undecapeptide ring.

Trifluoroacetic acid (TFA) is used as a reagent in the chromatographic purification of SCY-502635. The concentration of TFA in peptides is generally controlled at 1–10 mg/g [3], and the specification chosen for TFA in SCY-502635 is currently not more than 1 mg/g.

Previous studies concerning the assay of TFA in a variety of sample matrices have used headspace-GC [4], capillary electrophoresis [5] and, more commonly, ion chromatography [6–8]. Initially, a headspace-GC-FID method was developed and validated to control TFA in SCY-502635. This involved dissolving approximately 200 mg of drug substance with methanol, in the presence of sulphuric acid,

to convert the TFA present to the more volatile methyl ester (boiling points for TFA and its methyl ester are 72 °C and 43 °C, respectively; [4]). The limit of quantitation (LOQ) of this headspace-GC method was 0.05 mg/g. Although this degree of sensitivity was considered adequate for quality control testing (assuming a 1.0 mg/g specification), a lower limit of quantitation was desired for process monitoring and optimisation. It was also considered desirable, given the relatively high cost of synthesis and purification, to use a smaller sample mass for the TFA assay.

With these considerations in mind, ion chromatography was investigated as a possible alternative to GC. LOQ values reported in previous ion chromatography studies ranged from 0.4 ng to 10 ng TFA on-column. It is not possible to derive equivalent sensitivity data from the capillary electrophoresis study owing to the practical difficulty in measuring the injection volume, but the authors reported a concentration-based LOQ value of 1.2 µg/mL.

SCY-502635 is poorly water-soluble in acidic solution and practically insoluble at alkaline pH (see Table 1). It is, however, extremely soluble in methanol. Ion chromatography columns are prone to fouling by strongly-retained organic sample components [9,10], and it is important, in the interests of method robustness, either to limit the amount of these materials injected onto the column or to remove them from the column between sample injections using an organic eluent. The latter approach, however, has the disadvantage of extending the cycle time between injections owing to the need to re-equilibrate the column with aqueous mobile phase following the solvent rinse step. The Dionex AS11-HC ion

* Corresponding author. Tel.: +44 151 203 9812; fax: +44 1244 281283.

E-mail address: m.powell@quaypharma.com (M. Powell).

¹ Present address: Seachaid Pharmaceuticals, 801 Capitola Drive, Suite 5, Durham, NC 27713, United States.

² Present address: United Therapeutics Corporation, 55 T.W. Alexander Drive, Research Triangle Park, NC 27709, United States.

Table 1
Solubility of SCY-502635 at different pH values.

pH	Solubility (mg/mL)
1	103.66
3	40.59
5	12.03
7.5	1.2
8.6	0.62

exchange column was selected for this study owing to its high ion exchange capacity (290 μeq), which is ideal for the retention of singly-charged organic species, and its compatibility with organic solvents.

2. Materials and methods

2.1. Reagents

All reagents used were of analytical grade or better. Sodium trifluoroacetate of certified purity was supplied by Sigma–Aldrich Schweiz (Fluka Analytical), Buchs (St. Gallen), Switzerland. This reagent is hygroscopic and was stored in a desiccator when not in use. Carbonate-free sodium hydroxide (50% w/w aqueous solution) was purchased from the same supplier. Methanol (HPLC gradient grade) was supplied by Mallinckrodt Baker (Deventer, The Netherlands). Ultra-pure water of 18 M Ω resistivity or better was produced as required using a Purelab Option reverse-osmosis unit coupled to a Maxima water treatment system (Elga Process Water, Marlow, UK). Polyvinylidene fluoride (PVDF) membrane syringe filters (13 mm diameter, 0.2 μm porosity) were supplied by Whatman International (Maidstone, UK). A 19 mM solution of sodium hydroxide was prepared by dilution of a 50% (w/w) solution with ultra-pure water. Standards and spiking solutions were prepared by making suitable dilutions of an aqueous standard of sodium trifluoroacetate with ultra-pure water, the concentration of these solution being expressed in terms of the free acid.

2.2. Equipment

The ion chromatograph used for this study consisted of a GP50 quaternary gradient pump, AS50 autosampler, LC30 chromatography oven and an ED50 electrochemical detector (Dionex Limited, Camberley, UK). Eluents were degassed on-line using a series 200 vacuum degasser (Perkin-Elmer, Seer Green, UK) and hydroxide-containing mobile phase was maintained under nitrogen at approximately 3 psi to prevent ingress of atmospheric carbon dioxide. The autosampler was equipped with a 25 μL injection loop and the detector was fitted with a conductivity flow-cell. The background conductivity of the column effluent was reduced and the analyte signal enhanced using an ASRS 300 self-regenerating suppressor (Dionex, Camberley, UK) operated in recycle mode




Fig. 2. Overlay chromatogram showing a reagent blank (bottom), 0.1 $\mu\text{g/mL}$ TFA standard (middle) and 1.0 $\mu\text{g/mL}$ TFA standard (top). The TFA peak elutes at approximately 10.8 min.

using a current of 50 mA. The column used for the study was an IonPac AS11-HC hydroxide-selective anion exchange column (4 mm \times 250 mm; Dionex Part No. 52960) with a matching guard column (4 mm \times 50 mm, Part No. 52962).

2.3. Method development

2.3.1. Chromatographic conditions

Initial method conditions were assessed using 19 mM sodium hydroxide as the eluent at a flow rate of 1.0 mL/min. This eluent concentration was selected based on a previous study using the same type of ion chromatography column [6]. The column, conductivity flow-cell and suppressor were mounted in the chromatography oven and maintained at 30 $^{\circ}\text{C}$. The injection volume used was 25 μL , achieved by over-filling the injector loop. Sensitivity was evaluated by injecting standard solutions containing 0.1 $\mu\text{g/mL}$ and 1.0 $\mu\text{g/mL}$ TFA; these concentrations were close to the anticipated LOQ.

2.3.2. Sample preparation

The aim of sample preparation was to quantitatively recover TFA from the drug substance and to introduce the sample to the chromatographic system in such a way that robustness was not compromised. It was decided to try and avoid column clean-up between injections using an organic eluent since this would increase the time between sample injections. As has already been noted, SCY-502635 is virtually insoluble in alkaline aqueous solution (0.62 mg/mL) but readily soluble in methanol (>130 mg/mL). A sample preparation method was therefore evaluated based on complete dissolution of the sample in methanol followed immediately by dilution in mobile phase (19 mM sodium hydroxide) to remove SCY-502635 by precipitation. This approach also has the advantage of maximising the solubility and minimising the fugacity of TFA through the use of an alkaline sample diluent. Based on the instrument response from 0.1 $\mu\text{g/mL}$ and 1.0 $\mu\text{g/mL}$ chromatograms (Fig. 2), it was thought that dilution of a 50 mg sample of drug substance to a final volume of 25 mL would yield adequate sensitivity. Under these conditions, 1 mg/g TFA in the drug substance would produce a TFA concentration in the final extract of 2.0 $\mu\text{g/mL}$.

Four 50 mg aliquots of SCY-502635 were used to evaluate this sample preparation approach. Two of the samples were spiked with 50 μg TFA, equivalent to a concentration of 1 mg/g in the drug substance. All four samples were then dissolved in 1.0 mL methanol and 25.0 mL 19 mM sodium hydroxide was added. On addition of the hydroxide solution, a dense white precipitate was formed, which was removed by filtration through a 0.2 μm PVDF syringe filter. The potential for the filter to leach interfering ions




Fig. 1. SCY-502635.

Download English Version:

<https://daneshyari.com/en/article/1242348>

Download Persian Version:

<https://daneshyari.com/article/1242348>

[Daneshyari.com](https://daneshyari.com)