



Data in Brief

Genome-wide profiling to analyze the effects of Ox-LDL induced THP-1 macrophage-derived foam cells on gene expression



Yan-Wei Hu¹, Jia-Yi Zhao¹, Shu-Fen Li, Qian Wang, Lei Zheng*

Laboratory Medicine Center, Nanfang Hospital, Southern Medical University, Guangzhou, Guangdong 510515, China

ARTICLE INFO

Article history:

Received 22 September 2014

Accepted 29 September 2014

Available online 5 October 2014

Keywords:

Ox-LDL

Foam cell

LncRNA

Macrophages

Microarray

ABSTRACT

Atherosclerosis has a high incidence and is harmful to human health. An elevated level of oxidized low-density lipoprotein (Ox-LDL) is one of the major risk factors for atherosclerosis. During atherogenesis progression, circulating monocytes adhere to the intima and differentiate into macrophages. After differentiation, intimal macrophages intake Ox-LDL via scavenger receptors, thereby transforming into foam cells. Foam cell formation due to excessive accumulation of cholesterol by macrophages is a pathological hallmark of atherosclerosis. To gain a molecular understanding of the effect of Ox-LDL in atherosclerosis development, we conducted a genome-wide analysis of the Ox-LDL-induced macrophage transformation by microarray gene expression profiling. Here we describe in details the contents and quality controls for the gene expression and related results associated with the data uploaded to Gene Expression Omnibus (accession number GSE54039).

© 2014 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY-NC-SA license (<http://creativecommons.org/licenses/by-nc-sa/3.0/>).

Specifications

Organism/cell line/tissue	Homo sapiens
Strain(s)	Human monocytic THP-1 cells (Manassas, VA)
Sequencer or array type	Human LncRNA Array v2.0 (8 × 60 K, Arraystar)
Data format	Raw data: TXT files
Experimental factors	THP-1 macrophage cells without or with Ox-LDL treatment
Experimental features	Microarray gene expression profiling to identify genes that are regulated by Ox-LDL
Consent	N/A
Sample source	Location N/A

Direct link to deposited data

Deposited data are available here: <http://www.ncbi.nlm.nih.gov/geo/query/acc.cgi?acc=GSE54039>.

Experimental design, materials and methods

Cell line

The Human monocytic THP-1 cells were obtained from ATCC (Manassas, VA). It has been proved previously that THP-1 cells can

differentiate into macrophages following treatment with 100 nM phorbol 12-myristate 13-acetate (PMA) for 72 h and macrophages can transform into foam cells by incubation in the presence of 50 µg/ml oxidized low-density lipoprotein (Ox-LDL) for 48 h. Therefore, THP-1 macrophage-derived foam cells have been used as a model system for investigating the possible molecular mechanisms that Ox-LDL influences foam cell development, which represents a promising therapeutic strategy for atherosclerosis.

Ox-LDL was obtained from Biomedical Technologies Inc. (Stoughton, MA). Accumulated evidence showed that Ox-LDL played a significant role in the initiation and progression of atherosclerosis [1]. The previous study revealed that an elevated level of Ox-LDL was one of the major risk factors for atherosclerosis. Circulating monocytes could adhere to the intima and then differentiate into macrophages [2]. After differentiation, intimal macrophages intake Ox-LDL via scavenger receptors, thereby transforming into foam cells [3]. To gain a molecular understanding of the mechanism by which Ox-LDL transforms the macrophages into foam cells, we generated THP-1 macrophage-derived foam cells by Ox-LDL treatment.

Long noncoding RNAs (lncRNAs) are nonprotein coding transcribed RNA molecules consisting of more than 200 nucleotides [4]. An increasing number of studies have revealed that lncRNAs have a variety of important functions, including their role in transcription, splicing, translation, nuclear factor trafficking, imprinting, genome rearrangement, and chromatin modification [5,6]. Studies have shown that aberrant lncRNA expression is associated with a variety of human diseases such as cardiovascular diseases and cancer [7,8]. Thus, a better understanding of the roles of lncRNA in atherosclerosis will advance us to find a promising strategy for regulating lipid metabolism balance.

* Corresponding author. Tel./fax: +86 20 62787681.

E-mail address: nfyyl@163.com (L. Zheng).

¹ These authors contributed equally to this work.

Microarray and quality control

To identify the genes that are regulated by Ox-LDL, we isolated total RNA from 6 independent samples of THP-1 macrophages cultured in the presence or absence of 50 µg/ml of Ox-LDL in serum-free RPMI 1640 medium containing 0.3% BSA for 48 h using Trizol (Invitrogen Life Technologies). RNA from 6 samples was subjected to microarray analysis. Total RNA from each sample was quantified by the NanoDrop ND-1000 and RNA integrity was assessed by standard denaturing agarose gel electrophoresis (Fig. 1A and B). For microarray analysis, Agilent Array platform was employed. The sample preparation and microarray hybridization were performed based on the manufacturer's standard protocols with minor modifications. Briefly, mRNA was purified from 1 µg total RNA after removal of rRNA (mRNA-ONLY™ Eukaryotic mRNA Isolation Kit, Epicentre). Then, each sample was amplified and transcribed into fluorescent cRNA along the entire length of the transcripts without 3' bias utilizing a random priming method. The labeled cRNAs were hybridized onto the Human LncRNA Array v2.0 (8 × 60 K, Arraystar). After having washed the slides, the arrays were scanned by the Agilent Scanner G2505B. Agilent Feature Extraction software (version 10.7.3.1) was used to analyze acquired array images. All microarray data were submitted to the Gene Expression Omnibus (accession number GSE54039).

Normalization

Agilent Feature Extraction software (version 10.7.3.1) was used to analyze acquired array images. Quantile normalization and subsequent data processing were performed using the GeneSpring GX v11.5.1 software package (Agilent Technologies). After quantile normalization of the raw data, lncRNAs and mRNAs that at least 1 out of 6 samples have flags in Present or Marginal (All Target Values) were chosen for

further data analysis. Differentially expressed lncRNAs and mRNAs with statistical significance were identified through Volcano Plot filtering. Pathway analysis and GO analysis were applied to determine the roles of these differentially expressed mRNAs played in these biological pathways or GO terms. Finally, Hierarchical Clustering was performed to show the distinguishable lncRNA and mRNA expression pattern among samples.

Basic analysis

To explore possible changes in RNA expression during macrophage formation, we performed microarray analysis of THP-1 macrophages and THP-1 macrophage-derived foam cells using the Arraystar probe dataset, which included 24,748 lncRNAs and 24,420 coding transcripts. lncRNA and mRNA expression profiles from THP-1 macrophages (3 samples) and THP-1 macrophage-derived foam cells (3 samples) were produced using the Arraystar Human LncRNA Array v2.0 platform. Expression values were normalized based on the mean expression value for each probe set. Scatter plots of normalized signal values in lncRNAs and mRNAs are shown in Fig. 2A and B. Differently expressed probe sets were identified based on Student's *t*-test for paired samples' normalized expression values using the following cutoff: absolute fold change value larger than 3 and *P* value less than 0.01. By using a ≥3 fold change as a cut-off, *P*-value ≤0.01, 100 lncRNAs and 63 mRNAs were differentially expressed between THP-1 macrophage versus THP-1 macrophage-derived foam cells. Value distribution of signals from microarrays is shown in Fig. 3A and B.

Discussion

We described here a dataset composed of microarray expression profiling of Ox-LDL-responsive genes. With this dataset, we were able

A

Sample ID	OD260/280 Ratio	OD260/230 Ratio	Conc.(ng/µl)	Volume(µl)	Quantity(ng)	QC result Pass or Fail
Control 1	1.97	2.31	573.92	40	22956.80	pass
Control 2	1.98	2.32	751.90	30	22557.00	pass
Control 3	1.99	2.33	591.86	40	23674.40	pass
Treated 1	1.94	2.31	550.89	30	16526.70	pass
Treated 2	1.98	2.33	696.08	20	13921.60	pass
Treated 3	1.98	2.30	792.09	20	15841.80	pass



Lan1: Total RNA of sample Control 1
 Lan2: Total RNA of sample Control 2
 Lan3: Total RNA of sample Control 3
 Lan4: Total RNA of sample Treated 1
 Lan5: Total RNA of sample Treated 2
 Lan6: Total RNA of sample Treated 2

Fig. 1. Quality control of RNA. (A) RNA from 6 independent samples was quantified by NanoDrop ND-1000. (B) RNA integrity was assessed by standard denaturing agarose gel electrophoresis. Control 1 to control 3 represent RNA from THP-1 macrophages (three samples). Treated 1 to treated 3 represent RNA from THP-1 macrophage-derived foam cells (three samples).

Download English Version:

<https://daneshyari.com/en/article/2822291>

Download Persian Version:

<https://daneshyari.com/article/2822291>

[Daneshyari.com](https://daneshyari.com)