



Microbial Physiology

The dynamic structure of Spitzenkörper of *Trichosporon asahii* examined by the fluorescent probe FM4-64

Yuntong Lv[◇], Zhikuan Xia[◇], Wenling Wang*

Department of Dermatology, General Hospital of Beijing Military Command of PLA, 5 Nanmencang, Dongcheng District, Beijing 100700, PR China

ARTICLE INFO

Article history:

Received 11 April 2013

Accepted 17 June 2015

Associate Editor: Luis Henrique Souza Guimarães

Keywords:

Trichosporon asahii

Spitzenkörper

Polarized growth

Cytochalasin D

ABSTRACT

The Spitzenkörper is a dynamic and specialized multicomponent cell complex present in the tips of hyphal cells. The amphiphilic styryl dye FM4-64 was found to be ideal for imaging the dynamic changes of the apical vesicle cluster within growing hyphal tips. It is widely used as a marker of endocytosis and to visualize vacuolar membranes. Here we performed uptake experiments using FM4-64 to study the dynamic of the Spitzenkörper in *Trichosporon asahii*. We observed that Spitzenkörper were present at the tip of the budding site of the spore, blastospore, and the germ tube of *T. asahii*. We also found that Spitzenkörper were present at the tip of the hyphae as well as the subapical regions. Cytochalasin D, an inhibitor of actin polymerization, leads to abnormal Spitzenkörper formation and loss of cell polarity.

© 2015 Sociedade Brasileira de Microbiologia. Published by Elsevier Editora Ltda. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

Introduction

Studies have shown that vesicle trafficking is fundamental to numerous activities in eukaryotic organisms, and is involved in cell growth and differentiation. The vesicle trafficking network includes exocytosis and endocytosis.¹ In filamentous fungi, vesicle trafficking is concerned with the polarized tip growth.^{2,3} Vesicle trafficking to the apex during tip growth is highly organized and involves the activity of a specific, multicomponent organelle complex which, in most cases, is called the Spitzenkörper. The

Spitzenkörper is a vesicle-organizing center that is intimately involved in fungal cell morphogenesis, and it spatially and temporally coincides with polarity-related components. The Spitzenkörper was originally recognized as a dark region in phase-contrast microscopy at the tip of actively growing hyphae.⁴ Subsequently, electron microscopy revealed that this structure is rich in secretory vesicles.⁵ During tip growth, the dynamic behavior of the Spitzenkörper is predominated by secretory vesicles that form an “apical vesicle cluster,” which is intimately associated with the precise growth pattern of the hyphal apex.⁶

* Corresponding author.

E-mail: wangwl66@163.com (W. Wang).[◇] These authors contributed equally to this work.
<http://dx.doi.org/10.1016/j.bjm.2015.11.002>

The Spitzenkörper regulates hyphal apical growth primarily via exocytosis. The role of Spitzenkörper in endocytosis and recycling is currently not clear. Here we report that the Spitzenkörper comprises early compartments of the endocytic pathway in *Trichosporon asahii*. We performed uptake experiments using FM4-64, a widely used dye for endocytosis, to study the vesicle trafficking.⁷ We observed that cytochalasin D, an inhibitor of actin polymerization, leads to abnormal Spitzenkörper formation and loss of cell polarity.

Methods

Strains and growth conditions

The type strain of *T. asahii*, strain CBS 2479, was purchased from Centraalbureau voor Schimmelcultures (Utrecht, Netherlands) and used in this study. Cells were propagated by incubation in YPD (1% yeast extract, 2% peptone, 2% dextrose) medium at 37 °C. For FM4-64 (Biotium) staining, 5 µL of the stock solution (1.64 mmol/L in dimethyl sulfoxide (DMSO)) was added directly to 20 mL of culture, followed by incubation for an additional 10 min. The cells were harvested by brief centrifugation and washed with phosphate-buffered saline (PBS, pH 7.4).

Spitzenkörper inhibition

Cytochalasin D (Fermentek Ltd.) was dissolved in DMSO to obtain a concentration of 2 mmol/L, and this was added to culture to a final concentration of 0.2, 0.5, 1, 2, or 5 µmol/L. Same amount of DMSO without Cytochalasin D was used as a control. The effects of the various concentrations of inhibitor were tested by pretreating the samples with the indicated compound for 15 min prior to the addition of FM4-64. In some experiments, the cells were treated with the inhibitor for 2 h, after which the inhibitor solution was removed, and the cells were washed once and then resuspended in fresh medium without inhibitor or dye. The cells were then allowed to recover for up to 4 h, after which the dye was added to the cells for analysis.

Fluorescence microscopy

Mounted slides were photographed using a Leica DM 3000 fluorescence microscope with a Leica DFC300FX camera and an LP560 filter (Ex543/Em560). The viability of the single-celled propagules was confirmed by taking aliquots and determining the percentage that germinated after overnight incubation in the medium. Viability of hyphae was confirmed by their general appearance, presence of cytoplasmic streaming, and observation of continued tip elongation and growth. Images were recorded using a real-time digital imaging setup. At least 30 individual cells were examined for each experiment, and all the experiments were performed using at least three independent batches of cells.

Results

Uptake of FM4-64 at different stages of cell growth and development

FM4-64 is a lipophilic probe that fluoresces strongly after binding to the outer plasma membrane. Ten minutes after FM4-64 was added to the culture medium, fluorescence was observed in both apical and subapical hyphal compartments. To determine the kinetics of FM4-64 staining, cells were cultured at room temperature in the presence of FM4-64 and samples were examined at various time-points.

At 8 h, a single fluorescent spot with a diameter of 0.5–1 µm was observed at the tip of the budding site of spores (Fig. 1A). Similar fluorescent spots were observed at the tip of the germ tube (Fig. 1B) at 16 h. Dim background staining was observed in the cytoplasm (Fig. 1A and B). After 24 h of culture, a larger and brighter fluorescent spot (1–2 µm in diameter) was observed at the apical or subapical hyphae (Fig. 1C). This fluorescent spot resembles the Spitzenkörper structure.⁸ Fluorescent staining was also observed at septins (Fig. 1D), which could be the sites for branch projection. In addition, strong fluorescence was observed in the elongated spores of the pseudohyphae (Fig. 1E), as well as the chain-like spores (Fig. 1F). No Spitzenköpfer were noted in these structures.

Inhibition of Spitzenkörper formation

FM4-64 uptake was assessed in the presence of Cytochalasin D, an inhibitor of filamentous actin. *T. asahii* cells were treated with different concentration of Cytochalasin D (0.2, 0.5, 1, 2, or 5 µmol/L) for 15 min prior to the addition of the FM4-64, and the inhibitor was kept in the solution throughout the course of the experiment.

Light microscopy observation revealed that as the increase of Cytochalasin D concentration (from 0.2 to 5 µmol/L), the inhibition effects also increased. Hyphal growth became slower and hyphae became shorter but wider (Fig. 2).

When treated with a low dose of Cytochalasin D (0.2 µmol/L), the fluorescence was scattered in the central portion of the hyphae, but no obvious fluorescence was observed at the top of the hyphae (Fig. 3B). As the concentration increased, fluorescence was scattered further into the entire cell body (Fig. 3C–E). At 5 µmol/L, the hyphal growth was completely blocked and only a low background fluorescence was observed (Fig. 3F). No Spitzenkörper was observed when cells were treated with any concentration of Cytochalasin D. In the absence of Cytochalasin D, FM4-64-stained apical vesicle clusters at the tip of hyphae (Fig. 3A).

Discussion

The Spitzenkörper is a dynamic structure and closely related to polarized growth. It is continuously present at the hyphal tip at all stages of the cell cycle. In this study, we revealed the presence of Spitzenköpfer at the tip of the budding site of the spore, blastospore, and the germ tube of *T. asahii*. We also found that Spitzenköpfer were present at the tip of the

Download English Version:

<https://daneshyari.com/en/article/4356641>

Download Persian Version:

<https://daneshyari.com/article/4356641>

[Daneshyari.com](https://daneshyari.com)