

Contents lists available at ScienceDirect

Data in Brief

journal homepage: www.elsevier.com/locate/dib

Data Article

Primers and probe design and precision assessment of the real time RT-PCR assay in Coxsackievirus A10 and enterovirus detection



Jingfang Chen, Rusheng Zhang, Xinhua Ou, Dong Yao, Zheng Huang, Linzhi Li, Biancheng Sun*

Laboratory of Microbiology, Changsha Center for Disease Control and Prevention, Changsha, China

ARTICLE INFO

Article history: Received 7 February 2017 Received in revised form 7 April 2017 Accepted 25 April 2017 Available online 28 April 2017

Keywords: Real time RT-PCR Coxsackievirus A10 Enterovirus

ABSTRACT

This data article contains data related to the research article entitled "Rapid detection of enterovirus and Coxsackievirus A10 by a TaqMan based duplex one-step real time RT-PCR assay" (Chen at al., 2017) [1]. Primers and probe sequence design are among the most critical factors in real-time polymerase chain reaction (PCR) assay optimization. Linearity, sensitivity, specificity and precision are the crucial criteria which are used to evaluate the performance of a new method. This data article report the primers and probe design and precision assessment of the new assay. VP1 gene of Coxsackievirus A10 (CV-A10) and 5'-NCR of different enterovirus (EV) serotypes were retrieved from GenBank database and aligned. The intra- and inter-assay variation were assessed using high, medium and low concentration of control plasmid DNA and viral RNA samples.

© 2017 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY license (http://creativecommons.org/licenses/by/4.0/).

DOI of original article: http://dx.doi.org/10.1016/j.mcp.2017.02.003

* Corresponding author.

http://dx.doi.org/10.1016/j.dib.2017.04.035

2352-3409/© 2017 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY license (http://creativecommons.org/licenses/by/4.0/).

E-mail address: sunbiancheng2013@163.com (B. Sun).

Specifications Table

Subject area More specific	Biology Molecular Biology, real time RT-PCR
subject area	Molecular biology, fear time K1-PCK
Type of data	Table, figure
How data was	In silico analysis of gene sequences using online bioinformatics tools and
acquired	MEGA 5.2 software; Precision assay acquired by analysis of the threshold cycle
	value of control and clinical samples.
Data format	Raw, analyzed
Experimental	Gene sequences were retrieved from GenBank database; Standard plasmid
factors	DNA were constructed; Plasmid DNA and viral RNA concentration were
	quantified and the genome copies were calculated
Experimental	Primers and probe were designed using Primer Express software (version 3.0;
features	Applied Biosystems) and assay precision were determined by real time RT-PCR
Data source	Changsha, China
location	
Data accessibility	Data with this and the main article

Value of the data

- Rapid detection is crucial for Coxsackievirus A10 control along with the increasing circulation worldwide in the recent years.
- Specific and conserved regions for primers and probe design are defined using multiple sequence alignment.
- The intra- and inter-assay reproducibility are assessed using different concentrations of plasmid DNA and viral RNA.

1. Data

Using TaqMan probes we have previously established a real time RT-PCR method for Coxsackievirus A10 and other enterovirus detection [1].

The data presented in this article show the conserved regions for primers and probe design using multiple sequence alignment (Figs. 1 and 2). Table 1 represents data of the average threshold cycle (Ct) value, standard deviation (SD) and the coefficient of variability (CV) of different concentrations of plasmid DNA and viral RNA.

2. Experimental design, materials and methods

2.1. Primers and probes design

Sequences of CV-A10 VP1 gene and 5'-NCR of different enterovirus serotypes were downloaded from NCBI (http://www.ncbi.nlm.nih.gov/) database. Multiple sequences alignment was performed using the Clustal W algorithm in MEGA 5.2 software. Primer Express software (version 3.0; Applied Biosystems) was used to design the primers and probe for CV-A10 assay based on their highly conserved and specific regions (Fig. 1). For EVs detection, the primers and probe were chosen according to previous studies [2,3]; A small amplicon of 142 nucleotides in size were produced (Fig. 2). All oligonucleotides were synthesized by TAKARA (Dalian, China).

Download English Version:

https://daneshyari.com/en/article/4765180

Download Persian Version:

https://daneshyari.com/article/4765180

Daneshyari.com