



Copper-mediated N- and O-arylations with arylboronic acids in a continuous flow microreactor: a new avenue for efficient scalability

Brajendra K. Singh^{a,c}, Christian V. Stevens^{b,*}, Davy R. J. Acke^b, Virinder S. Parmar^c, Erik V. Van der Eycken^{a,*}

^aLaboratory for Organic and Microwave-Assisted Chemistry (LOMAC), Department of Chemistry, University of Leuven, Celestijnenlaan 200F, B-3001 Heverlee, Leuven, Belgium

^bResearch group SynBioC, Department of Organic Chemistry, Faculty of Bioscience Engineering, Ghent University, Coupure links 653, B-9000 Gent, Belgium

^cBioorganic Laboratory, Department of Chemistry, University of Delhi, Delhi 110 007, India

ARTICLE INFO

Article history:

Received 19 August 2008

Revised 23 September 2008

Accepted 26 September 2008

Available online 2 October 2008

Keywords:

Copper-mediated

N-Arylation

O-Arylation

Arylboronic acid

Continuous flow

Microreactor

Scalability

ABSTRACT

A continuous flow procedure has been elaborated for the copper(II)-mediated N- and O-arylation of a range of compounds with arylboronic acids using a commercial microreactor setup. The compounds could be continuously generated in good yields paving the way for efficient scalability.

© 2008 Published by Elsevier Ltd.

Among the C-heteroatom bond forming reactions, the copper-mediated heteroatom arylation has gained increasing importance. Independently reported by the groups of Chan,¹ Evans² and Lam³ in 1998, this cross-coupling procedure resolved the long-standing problem of the simple and mild formation of C(aryl)–O and C(aryl)–N bonds using an arylboronic acid as the aryl donor in combination with a base and copper acetate as the catalyst, as these reactions can be conducted at room temperature in air, providing a practical advantage over the Buchwald–Hartwig cross-coupling reaction.⁴ The method has proven to be applicable to a wide range of nucleophilic reaction partners including phenols, amines, anilines, amides, imides, ureas, carbamates, and sulfonamides. Also, aromatic heterocycles as imidazoles, pyrazoles, triazoles, tetrazoles, benzimidazoles, and indazoles participate efficiently in this copper(II)-mediated reaction.⁵ The scope of the procedure was expanded by the potential of a catalytic variant using less than stoichiometric amounts of Cu(OAc)₂ employing oxygen as oxidant or TEMPO/air, pyridinium N-oxide/air or [{Cu(μ-OH)(tmeda)}₂]Cl₂/oxygen in the presence of Cu(OAc)₂ as the Cu(II)-source.⁶

We recently reported the usefulness of this copper(II)-mediated cross-coupling protocol for the N-arylation of the 2(1H)-pyrazi-

none scaffold in solution phase⁷ and on solid support.⁸ As 2(1H)-pyrazinones have received considerable interest due to their valuable application as scaffolds for the generation of a diverse array of biologically interesting compound libraries,⁹ we were keen to know whether this cross-coupling protocol could be run in a continuous flow set up paving the way for efficient scalability. Microreactor technology is a rather new type of technology, which has gained high interest during the last two decades. The first literature related to chemistry in microstructured devices appeared in 1986 in a German patent that described how a microreactor should be built.¹⁰ Since then the interest in the application of microreactors is growing steadily.^{11,12} The main scale-up problems of organic reactions are associated with heat and mass transport, and shortcomings can lead to thermal runaway reactions. Microreactors should guarantee a better mass transfer, a safer handling of reactions, and a higher selectivity.¹³ However, the major advantage of this technology is the scalability.¹⁴ A classic scale-up is a time-consuming and expensive process. By contrast, once a process is optimized under microreactor conditions, there is no need anymore to scale up the reactor. Instead several microreactors are placed in parallel to produce higher quantities. Using this method, no modifications of the lab-scale conditions are necessary anymore.

Despite the large interest in copper-mediated heteroatom arylation during the last years,⁵ to our knowledge no studies are performed yet on the use of microreactors for this purpose. The microreactor that is used in our study is a commercially available

* Corresponding authors. Tel.: +32 16 327406.

E-mail addresses: chris.stevens@ugent.be (C. V. Stevens), erik.vandereycken@chem.kuleuven.be (E. V. Van der Eycken).




Figure 1. The applied microreactor: a CYTOS® College System, produced by CPC—Cellular Process Chemistry Systems GmbH; (1) input reagents via two pumps; (2) microreactor; (3) residence time unit; (4) output reaction mixture.

CYTOS® College System, produced by CPC—Cellular Process Chemistry Systems GmbH (Fig. 1).¹⁵ The entire device consists of two piston pumps, the microreactor itself and a residence time unit (RTU). The microreactor consists of several stacked plates. In a number of these plates, a path is etched. After stacking the plates together, a fluid pathway is formed with dimensions in the submillimeter range. The microreactor has two inlets. Each inlet flow is divided into four parallel flows which are one by one mixed together in an interdigital mixer. After passing through a mixing chamber, all parallel flows are once more divided into two flows each. Finally, all eight flows come together to leave the microreactor as one outlet flow (D).

Our initial investigations of this copper(II)-mediated cross-coupling in a microreactor setup were performed using pyrazinone **1** and phenylboronic acid (Scheme 1). For this purpose, a mixture of pyrazinone **1** (1 mmol, 238 mg) and boronic acid **2a** (2 mmol, 244 mg) in dichloromethane (DCM; 40 mL) was brought into the microreactor through one pump and a solution of Cu(OAc)₂ (1 mmol, 182 mg), Et₃N (1 mmol, 102 mg), and pyridine (2 mmol, 160 mg) in DCM (40 mL) through another pump allowing efficient mixing at rt in the microcapillary (Fig. 1). Different flow rates resulting in different residence times ranging from 15 to 120 min were applied (Table 1). When the reaction was performed using a relatively high flow rate of 1.54 mL/min, approximately 30% conversion to the desired N-arylated pyrazinone was obtained after 15 min (Table 1, entry 1). Prolonging the residence time to 45 min by lowering the flow rate to 0.52 mL/min resulted in an increased conversion of 50% (Table 1, entry 2). Consequently, when the reaction was carried out for 120 min, a 100% conversion was achieved delivering 79% of the desired N-arylated compound **3** after purification (Table 1, entry 5). This is much faster compared to the required 12 h under conventional conditions.⁷ Alternative solvents were screened, but only DMF seemed to be as efficient as dichloroethane (DCE); acetonitrile and dioxane were unsuccessful to dissolve the copper acetate efficiently resulting in clogging

Table 1

Optimization of the flow rate for the N-arylation of pyrazinone **1** using microreactor technology

Entry ^a	Residence time (min)	Flow rate (mL/min)	Conversion ^b (%)
1	15	1.54	30
2	45	0.52	50
3	60	0.39	70
4	90	0.26	90
5	120	0.20	100 ^c


^a All optimization reactions were performed at 25 °C on a 1.0 mmol scale of pyrazinone **1**. The pyrazinone **1** and boronic acid **2a** (2 equiv) were dissolved in DCM (40 mL) and Et₃N (1.0 equiv), pyridine (2.0 equiv) and Cu(OAc)₂ (1.0 equiv) were separately dissolved in DCM (40 mL), and both mixtures were simultaneously brought into the microcapillary by the aid of two separate pumps (with the same pump rate).

^b The conversion was estimated by quantitative TLC.

^c This resulted in an isolated yield of 79%.

problems of the microreactor. In DMF, the reaction was faster as compared to DCM as the reaction could be run at elevated temperature. When the reaction was performed at 130 °C for 30 min, a 100% conversion was achieved delivering compound **3** in 72% yield. However, isolation from the DMF mixture turned out to be difficult. Therefore, we decided to perform our further investigations using DCM at rt. We then gradually increased the concentration of different compounds to 5- and 10-fold, applying the optimized conditions, that is, a flow rate of 0.2 mg/mL in DCM at rt. Both reactions were running well without any problem. Applying the 10-fold concentration, that is, reaction on 10.0 mmol of pyrazinone **1**, we were able to continuously generate compound **3** in 21 mg (0.066 mmol) per minute, opening the way for its synthesis on a multigram scale. However, further increasing the concentration to a 15-fold (15 mmol of pyrazinone **1**) caused clogging problems of the microreactor due to the low solubility of Cu(OAc)₂ in DCM. To the best of our knowledge, copper(II)-mediated C–N coupling reactions have never been performed on a multigram scale applying microreactor technology.

After having established an efficient continuous flow microreactor protocol¹⁶ for the N-arylation of pyrazinone **1**, we evaluated the applicability for differently substituted anilines (Scheme 2, Table 2). Thus, aniline **4a** (1.0 mmol, 93 mg) and phenylboronic acid **2a** (2.0 mmol, 244 mg) were dissolved in DCM (40 mL), and mixed in the microcapillary with a solution of Cu(OAc)₂ (1.0 mmol, 182 mg), Et₃N (1.0 mmol, 101 mg), and pyridine (2.0 mmol, 158 mg) in DCM (40 mL) applying both pumps with a flow rate of 0.2 mL/min and a residence time of 2 h. The reaction proceeded well yielding 71% of the N-arylated aniline **5a** (Table 2, entry 1). The concentration could easily be increased to 10-fold resulting in the formation of **5a** at a rate of 10 mg/min. As a proof of concept, when the reaction was run under conventional conditions applying the same parameters, the arylation of 1 g of aniline with phenylboronic acid yielded the N-arylated aniline in 67% in 24 h clearly indicating the usefulness of our microreactor protocol. However, upon further scale up to 10-fold, we noticed that the reaction could not



Scheme 1. N-Arylation of the pyrazinone **1** scaffold.

Download English Version:

<https://daneshyari.com/en/article/5276257>

Download Persian Version:

<https://daneshyari.com/article/5276257>

[Daneshyari.com](https://daneshyari.com)