### Accepted Manuscript

Title: Examination of dual-pendant sulfonate bis(indolenine)squaraine with 1,8-diazabicyclo[5.4.0]undec-7-enium cations



Authors: Daniel E. Lynch, Martin Cox, Patrick M. King

 PII:
 S1010-6030(17)30391-X

 DOI:
 https://doi.org/10.1016/j.jphotochem.2017.06.032

 Reference:
 JPC 10711

To appear in: Journal of Photochemistry and Photobiology A: Chemistry

 Received date:
 22-3-2017

 Revised date:
 12-6-2017

 Accepted date:
 22-6-2017

Please cite this article as: Daniel E.Lynch, Martin Cox, Patrick M.King, Examination of dual-pendant sulfonate bis(indolenine)squaraine with 1,8-diazabicyclo[5.4.0]undec-7-enium cations, Journal of Photochemistry and Photobiology A: Chemistry https://doi.org/10.1016/j.jphotochem.2017.06.032

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.

## ACCEPTED MANUSCRIPT

## Examination of dual-pendant sulfonate bis(indolenine)squaraine with 1,8diazabicyclo[5.4.0]undec-7-enium cations

Daniel E. Lynch <sup>a,\*</sup>, Martin Cox <sup>b</sup> and Patrick M. King <sup>b</sup>

<sup>a</sup> Exilica Limited, The Technocentre, Puma Way, Coventry CV1 2TT, UK

<sup>b</sup> Faculty of Health and Life Sciences, Coventry University, Coventry CV1 5FB, UK

\* Corresponding author. Tel.: +44(0)24 7765 8505. Email address: d.lynch@exilica.co.uk

(D.E. Lynch).

#### GRAPHICAL ABSTRACT



#### HIGHLIGHTS

- 1,8-diazabicyclo[5.4.0]undec-7-ene used to make bis(indolenine)squaraine dyes
- Solution protein fluorescence responses comparable to bis-piperdinium analogue
- X-ray crystal structure of bis-(1,8-diazabicyclo[5.4.0]undec-7-enium) dye reported
- Dye-stained electrophoresis gels do not fluoresce
- Bacterial cells stained with bis-piperidinium compound fluoresce.

#### ABSTRACT

The dual-pendant sulfonate bis(indolenine)squaraines with 1,8diazabicyclo[5.4.0]undec-7-enium cations has been prepared and examined with reference to the previously studied bis-piperidinium analogue. The X-ray crystal structure of the bis-(1,8Download English Version:

# https://daneshyari.com/en/article/6492511

Download Persian Version:

https://daneshyari.com/article/6492511

Daneshyari.com