

Accepted Manuscript

Title: A novel method for synchronizing motion capture with other data sources for millisecond-level precision

Author: V. Komisar A.C. Novak Dr. B. Haycock

PII: S0966-6362(16)30603-8

DOI: <http://dx.doi.org/doi:10.1016/j.gaitpost.2016.10.002>

Reference: GAIPOS 5181

To appear in: *Gait & Posture*

Received date: 14-1-2016

Revised date: 26-9-2016

Accepted date: 2-10-2016



Please cite this article as: Komisar V, Novak AC, Haycock B. A novel method for synchronizing motion capture with other data sources for millisecond-level precision. *Gait and Posture* <http://dx.doi.org/10.1016/j.gaitpost.2016.10.002>

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.

A novel method for synchronizing motion capture with other data sources for millisecond-level precision

Komisar V^{1,2}, Novak AC^{1,3,4}, Haycock B¹

¹ iDAPT Centre for Rehabilitation Research, Toronto Rehabilitation Institute – University Health Network, 550 University Ave, Toronto, ON M5G 2A2, Canada, www.idapt.com

² Institute of Biomaterials and Biomedical Engineering, University of Toronto, 164 College Street, Toronto, ON, M5S 3G9, Canada

³ Department of Occupational Sciences and Occupational Therapy, University of Toronto, 500 University Ave, Toronto, ON, M5G 1V7, Canada

⁴ Faculty of Kinesiology and Physical Education, University of Toronto, 55 Harbord St, Toronto, ON, M5S 2W6, Canada

Corresponding Author:

Dr. Bruce Haycock

iDAPT (www.idapt.com), Toronto Rehabilitation Institute – University Health Network

550 University Ave, Toronto, Ontario

M5G 2A2 Canada

Email: bruce.haycock@uhn.ca

Tel: 416-597-3422 x7866

Fax: 416-597-3027

Highlights

- A simple method to align motion capture with other data sources is described
- A common analog signal enables continuous synchronization throughout entire trials
- Millisecond-level precision enables accurate interpretation of rapid movements
- Data from individual trials can be synchronized automatically

ABSTRACT

Download English Version:

<https://daneshyari.com/en/article/8798711>

Download Persian Version:

<https://daneshyari.com/article/8798711>

[Daneshyari.com](https://daneshyari.com)