Accepted Manuscript

Title: Direct detection of Tritrichomonas foetus in cattle genital fluid trough loop mediated isothermal amplification of elongation factor 1 alpha 1

Author: Jorge Oyhenart

PII: S0304-4017(18)30300-5

DOI: https://doi.org/10.1016/j.vetpar.2018.08.011

Reference: VETPAR 8736

To appear in: Veterinary Parasitology

Received date: 22-6-2018 Revised date: 6-8-2018 Accepted date: 24-8-2018

Please cite this article as: Oyhenart J, Direct detection of Tritrichomonas foetus in cattle genital fluid trough loop mediated isothermal amplification of elongation factor 1 alpha 1, *Veterinary Parasitology* (2018), https://doi.org/10.1016/j.vetpar.2018.08.011

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.



ACCEPTED MANUSCRIPT

Direct detection of Tritrichomonas foetus in cattle genital fluid trough loop mediated isothermal amplification of elongation factor 1 alpha 1

Jorge Oyhenart*

INCITAP - CONICET – Facultad de Ciencias Exactas y Naturales, Universidad Nacional de La Pampa, Av. Uruguay 151 (6300) Santa Rosa, La Pampa, Argentina.

*corresponding author: e-mail: jorgeoyhenart@gmail.com

Highlights

- Fast and sensitive method for *Tritrichomonas foetus* detection.
- Based on loop mediated isothermal amplification of *T. foetus* elongation factor 1 alpha 1.
- 100 times more sensitive than culture, tf-5.8r LAMP or TFR3-TFR4 PCR
- Direct detection from smegma and cervical mucus.

Download English Version:

https://daneshyari.com/en/article/9954686

Download Persian Version:

https://daneshyari.com/article/9954686

Daneshyari.com